



Cite this: DOI: 10.1039/d3sm01641b

## Enzymatic cleaving of entangled DNA rings drives scale-dependent rheological trajectories†

 Philip Neill,  Natalie Crist, Ryan McGorty  and Rae Robertson-Anderson \*

DNA, which naturally occurs in linear, ring, and supercoiled topologies, frequently undergoes enzyme-driven topological conversion and fragmentation *in vivo*, enabling it to perform a variety of functions within the cell. *In vitro*, highly concentrated DNA polymers form entanglements that yield viscoelastic properties dependent on the topologies and lengths of the DNA. Enzyme-driven alterations of DNA size and shape therefore offer a means of designing active materials with programmable viscoelastic properties. Here, we incorporate multi-site restriction endonucleases into dense DNA solutions to linearize and fragment circular DNA molecules. We pair optical tweezers microrheology with differential dynamic microscopy and single-molecule tracking to measure the linear and nonlinear viscoelastic response and transport properties of entangled DNA solutions over a wide range of spatiotemporal scales throughout the course of enzymatic digestion. We show that, at short timescales, relative to the relaxation timescales of the polymers, digestion of these ‘topologically-active’ fluids initially causes an increase in elasticity and relaxation times followed by a gradual decrease. Conversely, for long timescales, linear viscoelastic moduli exhibit signatures of increasing elasticity. DNA diffusion, likewise, becomes increasingly slowed, in direct opposition to the short-time behavior. We hypothesize that this scale-dependent rheology arises from the population of small DNA fragments, which increases as digestion proceeds, driving self-association of larger fragments *via* depletion interactions, giving rise to slow relaxation modes of clusters of entangled chains, interspersed among shorter unentangled fragments. While these slow modes likely dominate at long times, they are presumably frozen out in the short-time limit, which instead probes the faster relaxation modes of the unentangled population.

 Received 5th December 2023,  
 Accepted 27th February 2024

DOI: 10.1039/d3sm01641b

[rsc.li/soft-matter-journal](https://rsc.li/soft-matter-journal)

## Introduction

The cellular environment exists in a state far from equilibrium – a necessity to perform complex cellular functions.<sup>1</sup> Densely crowded biopolymers such as proteins and DNA are constantly manipulated by molecular motors and enzymes which drive changes to the architecture and transport of the biopolymers, which in turn alter the viscoelastic properties of the cellular environment.<sup>2–6</sup>

Much work over the last two decades has been devoted to how the action of molecular motors (*e.g.*, myosin, kinesin) on cytoskeleton filaments (*e.g.*, actin, microtubules) drive mechanical changes to *in vitro* cytoskeletal networks, as well as the cellular environment.<sup>7–11</sup> In these systems, the motors consume energy to actively pull on the filaments to drive mechanical changes.

However, only recently has the action of restriction enzymes, that alter the topology and/or length of biopolymers such as

DNA, been appreciated as a potential driver of changes to viscoelastic properties.<sup>12–18</sup> Yet, this process is ubiquitous in diverse processes including nuclear division, DNA repair, transcription and signaling.<sup>1,3,19,20</sup> For example, DNA naturally occurs in supercoiled circular, relaxed circular (ring), and linear topologies with lengths that scale many decades. A host of restriction endonucleases naturally exist to cleave DNA strands at specific sites, *via* basepair recognition, to convert supercoiled or ring DNA to linear form, with a single cut, and subsequently cleave the linear strand into smaller fragments if multiple recognition sites exist along the chain.<sup>19</sup>

Recent work demonstrated that topology-altering enzymatic digestion of DNA can push concentrated DNA solutions out of equilibrium to elicit time-varying microrheological properties of the solutions en route from their initial, undigested steady-state to their fully-digested steady-state.<sup>13</sup> In these studies, the viscosity of ‘topologically-active’ solutions of 5.9 kilobasepair (kbp) DNA was shown to steadily increase under the action of a single-site restriction enzyme that cleaves each circular DNA molecule into a linear strand of the same length. The degree to which the viscosity increased during linearization depended on the enzyme:DNA stoichiometry as well as the concentration of

Department of Physics and Biophysics, University of San Diego, San Diego, CA 92104, USA. E-mail: [randerson@sandiego.edu](mailto:randerson@sandiego.edu)

† Electronic supplementary information (ESI) available. See DOI: <https://doi.org/10.1039/d3sm01641b>

the DNA. This effect was shown to arise from the changing degree of polymer overlap as circular DNA (rings and supercoiled constructs) are converted to linear chains. Namely, circular polymers have smaller radius of gyration compared to their linear counterparts, following the respective empirical relations  $R_{G,L} \simeq 1.58R_{G,R}$  and  $R_{G,L} \simeq 2.1R_{G,S}$  for ring and supercoiled circular DNA.<sup>21</sup> Because the polymer overlap concentration  $c^* \simeq \left[ \frac{M}{N_A} \right] / \left[ \frac{4}{3} \pi R_G^3 \right]$  scales as  $R_G^{-3}$ ,<sup>22</sup> as circular DNA is converted to linear topology, the degree of polymer overlap  $c/c^*$  increases despite the fact that  $c$  does not change.

This same study showed that fragmentation of entangled linear DNA caused the viscosity to steadily decrease, as expected as the degree of polymer overlap also decreases with decreasing polymer length  $L$ .<sup>13</sup> Specifically, given the molecular weight of a polymer scales with its length  $M \sim L$  and for an ideal random coil  $R_G \sim L^{1/2}$  it follows that  $c^* \sim L^{-1/2}$  such that  $c/c^* \sim L^{1/2}$ .<sup>23</sup> The dependence of overlap and entanglements on length can also be seen in the reptation model relation between polymer length  $L$  and the polymer length between entanglements  $l_e \sim c^{-1}$ .<sup>31,22</sup> Because  $l_e$  is set by the concentration, the shorter the chains are at a given concentration the fewer entanglements per chain ( $L/l_e$ ), which dictates the viscoelastic response for entangled polymers.

Following on this work, bulk rheology studies of topologically-active 11 kbp ring DNA undergoing digestion by a 13-site cutter showed an unexpected sharp transition from a low-elasticity state to a high-elasticity state, marked by a  $\sim 50$ -fold increase in the elastic modulus  $G'$ , upon complete digestion.<sup>14</sup> Not only is the rise in elasticity at odds with the decreasing viscosity measured for fragmenting linear chains described above, but the transition between the states was also not steady or continuous, as was reported for microrheology measurements.<sup>13</sup> This same study also showed that the same enzymatic activity elicited very different rheological trajectories for composites of DNA and dextran with varying volume fractions of DNA. Namely, while all systems displayed the initial abrupt jump from low-elasticity to high-elasticity, this transition occurred much earlier in the digestion process for DNA-dextran composites, as compared to 100% DNA, and was followed by a subsequent drop back to a low-elasticity state. These results were shown to arise from cooperative percolation of clusters of slow (long, entangled) DNA, facilitated *via* crowding by small dextran polymers, that dictate the rheological response.<sup>14,24–26</sup>

The studies described above highlight the different responses that topologically-active DNA systems can exhibit at different spatiotemporal scales; and the important role that microscale heterogeneities, crowding, and depletion interactions may play. Moreover, the effect of end closure of ring polymers on their ability to form entanglements and elicit elastic-like response, a topic of intense interest and debate in the polymer community,<sup>24,27–41</sup> also contributes non-trivially to the rheological trajectory of these topologically-active systems. Several theoretical frameworks have been developed to model the steady-state steric interactions between highly overlapping


rings,<sup>27,33,34,38,41–45</sup> akin to entanglements for linear chains, treating rings as branched, amoeba-like, tree-like or loopy structures that exhibit weaker constraints and decreased elastic contribution to the stress response as compared to linear chains.<sup>31,46–48</sup> However, threading of rings by linear chains in transient ring-linear blends that form during digestion may exhibit even stronger constraints and slowed relaxation compared to transient 100% linear solutions, as evidenced in numerous theoretical and experimental investigations of steady-state ring-linear blends.<sup>30,34–36,40,47–55</sup> Moreover, recent simulations have shown evidence of ring–ring threadings in highly overlapping ring polymer melts, which leads to heterogeneous dynamics and cooperative glassiness,<sup>24,26,51,56</sup> consistent with macrorheology experiments.<sup>14</sup>

While the microscale and macroscale responses of topologically-active fluids are clearly different, the elusive mesoscale rheology that connects these two limits remains an important unknown. How the diffusion of the molecules undergoing digestion varies to give rise to the non-equilibrium rheological response, not yet understood, is critical to mapping the dynamics across spatiotemporal scales. Moreover, how the time-varying topologies of the molecules effect their response to nonlinear straining, at rates much faster than the intrinsic relaxation rates of the system, is also an essential unanswered question with a likely non-trivial answer.

For example, nonlinear optical tweezers microrheology experiments on blends of ring and supercoiled DNA have shown evidence of strain-induced separation of rings and supercoiled constructs, with rings building up at the leading edge of the trapped bead as it is dragged through the solution, while supercoiled constructs disentangle from the rings and fall into the wake of the strain path.<sup>57</sup> A similar phenomenon was observed in ring-linear DNA blends in which entangled linear chains built-up at the leading edge, increasing the local entanglement density, while rings fell into the wake.<sup>48</sup> Finally, simulations of entangled rings have demonstrated strain-induced deadlocking and enhanced threading of rings, forming superstructures and slow-moving clusters that lead to thickening and glassy dynamics not seen with entangled linear polymers.<sup>37,39</sup>

Here, we design topologically-active solutions of ring and supercoiled plasmid DNA of length  $L = 5.9$  kbp at a concentration well above the overlap concentration of the initially circular DNA,  $c/c^* \simeq 8$ .<sup>21</sup> As the DNA becomes linearized during digestion, the degree of overlap increases, approaching an upper bound of  $c/c_L^* \simeq 38$  for the case of a solution of fully linearized but not yet fragmented molecules.<sup>21</sup> This concentration equates to  $c/c_e \simeq 6.3$ , where  $c_e \simeq 6c_L^*$  is the critical entanglement concentration, and an entanglement density of  $Z \simeq (c/c_e)^{5/4} \simeq 10$ .<sup>23,58,59</sup> We incorporate a multi-site restriction endonuclease HaeIII, that cleaves the DNA at 17 different recognition sites, creating 17 blunt-ended linear strands ranging in length from 19 bp to 1642 bp (Fig. 1A and S1, ESI<sup>†</sup>).

To address the issues raised above, we measure the time-resolved rheological properties of the DNA solutions across three decades of spatiotemporal scales from  $t \simeq 17$  ms to



**Fig. 1** Optical tweezers microrheology and differential dynamic microscopy (DDM) couple time-dependent rheology with transport properties of topologically-active DNA across spatiotemporal scales. (A) Digestion of 5.9 kbp plasmid DNA (pYES2) by restriction enzyme HaeIII first linearizes the initially circular DNA molecules, then successively cleaves the full-length linear chains into shorter fragments over activity time  $t_a$ , resulting in 17 linear fragments of lengths  $l_i = 19$  bp to 1642 bp at the end of digestion. (B) and (C) Cartoon of optical tweezers microrheology (OTM) probing both the linear (B) and nonlinear (C) force response of concentrated DNA solutions as they undergo HaeIII-driven topological digestion. (B) A  $4.5 \mu\text{m}$  microsphere probe (grey) embedded in a DNA solution (green) and trapped by an infrared laser (red). The thermal fluctuations of the trapped probe from the surrounding DNA solution are measured by recording the corresponding deflections of the trapping laser via a position sensing detector (shown in inset in (D)). (C) The same microsphere is dragged through the DNA solution at a constant speed of  $v = 60 \mu\text{m s}^{-1}$  (left) through a distance  $s = 30 \mu\text{m}$ , using a nanopositioning stage to move the sample relative to the trap. After this displacement, the stage motion is halted, and relaxation of the force on the probe from the surrounding DNA relaxing strain-induced stress is measured (right). (D) Thermal fluctuations measured at 20 kHz over the course of 100 s (shown in inset) are used in linear OTM to compute the frequency-dependent linear viscoelastic moduli,  $G'(\omega)$  (filled squares) and  $G''(\omega)$  (open squares), using the generalized Stokes-Einstein relation (GSER), as described in Methods. (E) and (F) In nonlinear OTM, the force exerted on the probe by the DNA is measured during (E) and following (F) the  $30 \mu\text{m}$  displacement to determine the strain response and relaxation dynamics. (G) To measure DNA diffusion using DDM, a small fraction of DNA is fluorescently-labeled and time-series of the diffusing molecules are recorded. (H) DDM analysis results in an image structure function  $D(\vec{q}, \Delta t)$  as a function of lag time  $t$  and wave vector  $q$  in Fourier space. Shown is  $D(\vec{q}, \Delta t = 20 \text{ s})$ . The radially symmetric image structure function is azimuthally averaged to arrive at a one-dimensional function  $D(q, t)$  from which we extract the intermediate scattering function (ISF). (I) Sample ISF as a function of lag time  $t$  for  $q = 2 \mu\text{m}^{-1}$ , which is fit to an exponential to determine the (J) density fluctuation decay time  $\tau$  which exhibits power-law scaling  $\tau(q) \sim q^{-2}$ , indicative of diffusive dynamics.

$t \simeq 6.3 \text{ s}$  and from  $x \simeq 20 \text{ nm}$  to  $x \simeq 30 \mu\text{m}$ . We achieve this breadth by combining both linear and nonlinear optical tweezers microrheology measurements (OTM) (Fig. 1B–F) with differential dynamic microscopy (DDM) analysis of DNA fluctuations (Fig. 1G–J).

We perform measurements for two different HaeIII:DNA stoichiometries of  $0.3 \text{ U } \mu\text{g}^{-1}$  and  $0.5 \text{ U } \mu\text{g}^{-1}$ , chosen based on previous results, to be high enough to achieve complete digestion over the course of  $< 4$  hours but low enough such that the system can be considered in quasi-steady state during each

individual measurement. Specifically, the longest measurement timescale we probe ( $<10$  s) is several orders of magnitude shorter than the time to completely digest the DNA,  $t_d \approx 6000$  s, and time for a single cut  $t_{d,1} \approx t_d/17 \approx 350$  s. This separation of timescales, between measurement times  $t$  and enzymatic activity time  $t_a = 10$ –240 min, allows us to treat the solutions during each measurement as in quasi-steady state, such that we can evaluate steady-state ensemble-averaged quantities such as viscoelastic moduli,  $G'(\omega)$  and  $G''(\omega)$ , and diffusion coefficients  $D$ .

Throughout the paper, we will refer to the time associated with digestion as activity time  $t_a$ , the variation in measured quantities as a function of  $t_a$  as trajectories, and the timescales probed in each measurement as the measurement timescale or time  $t$ . While we acknowledge that enzymatic digestion is a stochastic process, so there may be topological and conformational heterogeneities that occur on timescales below  $t_{d,1}$ , we still expect our adiabatic assumption to be valid for understanding the ensemble dynamics we present.

Our results uncover stark differences between the trajectories of rheological metrics that probe short *versus* long measurement timescales. Namely, the response of topologically-active solutions to fast strain rates (short times) appear to first increase in elastic-like contributions followed by a large drop before reaching steady-state force and viscosity values that are substantially lower than the solution before digestion. In direct opposition, the slow strain rate (long-time) response and diffusivity initially show reduced elastic-like contributions and constraints to motion, followed by a surprising increase to values that are significantly higher (*i.e.*, more elastic-like) than their pre-digestion counterparts.

Our results shed important light on the interplay between topology and rheology at different scales, and highlight the unique heterogeneities and hierarchies that arise in topologically-active solutions due to time-varying contributions of topology, crowding, entanglements, threadings, and polymer overlap. These out-of-equilibrium fluids may have applications in timed self-assembly, selective filtration and sequestration, and wound-healing.

## Results

### Enzymatic digestion of circular DNA initially increases elastic contributions to the nonlinear force response followed by a decrease to a more viscous-dominated regime

We perform nonlinear optical tweezers microrheology (OTM) on enzymatically active solutions by pulling a bead of radius  $R = 2.25$   $\mu\text{m}$  a distance  $s = 30$   $\mu\text{m}$  through the solution at a constant speed of  $v = 60$   $\mu\text{m s}^{-1}$ . Following this strain, we hold the bead fixed and measure the force decay as the system relaxes into a new mechanically-steady state. The distance and speed were chosen to both be well above the intrinsic lengthscales and relaxation rates of the system to probe the nonlinear response. Namely, the fully extended length of the plasmid is  $L \approx 1.97$   $\mu\text{m}$  and the radius of gyration of the linear, ring, and supercoiled topologies are  $R_{G,L} \approx 213$  nm,  $R_{G,R} \approx 135$  nm, and  $R_{G,S} \approx 103$  nm.<sup>21</sup> Moreover, the strain speed can be approximately converted to a strain rate of  $\dot{\gamma} \approx \frac{3v}{\sqrt{2}R} \approx 57$   $\text{s}^{-1}$ ,<sup>60</sup>

which is faster than the intrinsic Rouse relaxation rates of  $\tau_{R,L}^{-1} \approx 7.4$   $\text{s}^{-1}$ ,  $\tau_{R,R}^{-1} \approx 24$   $\text{s}^{-1}$  and  $\tau_{R,S}^{-1} \approx 52$   $\text{s}^{-1}$  for, respectively, linear, ring, and supercoiled chains, computed from the expression  $\tau_R^{-1} \approx \pi D_0/2NR_G^2$  where  $N = 19.6$  and  $D_0$  is the topology-dependent tracer diffusion coefficient.<sup>21–23,61</sup>

As shown in Fig. 2A, the force response to nonlinear straining at all times exhibits an initial increase in force followed by a rolling over to a strain-independent plateau, indicative of a viscous response. The average value of this plateau force  $F_p$  increases in the first 10 min after digestion begins ( $t_a = 10$  min), followed by large drops from  $t_a = 20$  min to  $t_a = 50$  min, then a more modest decrease before reaching a  $t_a$ -independent plateau starting at  $t_a \approx 100$  min (Fig. 2B). Moreover, the force response for the initial time point ( $t_a = 10$  min) exhibits a stress overshoot, seen in macrorheology measurements of highly entangled linear polymers,<sup>22,62</sup> but typically not seen in microrheology measurements,<sup>63,64</sup> demonstrating strong steric constraints at this early stage of the digestion process. We observe similar, albeit weaker, overshoots at  $t_a \approx 0$  and at  $t_a \approx 20$  min, indicating that the underlying mechanism driving the overshoot is present in the initial purely circular DNA solution, but is enhanced by the introduction of linear chains into the circular solution at early digestion times, and then suppressed as digestion proceeds and circular DNA is eliminated. This trend suggests that threading events, which can occur between neighboring rings but are weaker than ring-linear threadings, may be the driver of the overshoot we observe which is eliminated at an activity time that we anticipate correlates with complete linearization and beginning of fragmentation.

The general trend of initially increasing in force then decreasing aligns with the results of ref. 65 if we attribute the initial rise to the linearization of rings, which increases polymer overlap and entanglements, followed by the drop that was previously seen for viscosity of linear chains undergoing fragmentation. In these previous works, evidence of ring-linear threading was also reported during the linearization of rings in solutions of the same concentration and length we study here. Similarly, and given the potential indication of threading described above, we also expect that ring-linear threading contributes to the initial rise in force and that the subsequent reduction occurs both from fragmentation but also by eliminating threading events.

Moreover, the functional form of the force decay over time and the final value of the plateau force is relatively insensitive to enzyme:DNA stoichiometry, with the lower value ( $0.3$   $\mu\text{g}^{-1}$ ) reaching only a slightly higher initial plateau force  $F_p(t_a = 10$  min) and approaching the terminal  $F_p$  value slightly later in time ( $\sim 100$  min *vs.*  $\sim 75$  min). This dependence aligns with a modestly slower digestion rate, which is predicted to scale with stoichiometry, but with the same mechanisms driving the rheological changes.

We also examine the small-strain behavior, during the initial elastic-like increase in force ( $0 < x < 1$   $\mu\text{m}$ ), to determine how the elastic response varies during digestion (Fig. 2C). As digestion proceeds ( $t_a$  increases), the force curves become shallower, indicating reduced elastic-like contributions; and




Fig. 2 Nonlinear OTM reveals transient increase followed by gradual decrease of resistive force and relaxation timescales of topologically-active circular DNA undergoing linearization and fragmentation. (A) Force exerted by the DNA as the trapped probe is moved relative to the sample, by increasing the  $x$ -position of the stage, for increasing activity times during HaeIII digestion.  $t_a = 0$  denotes the time at which  $0.5 \text{ U } \mu\text{g}^{-1}$  of HaeIII was added to the DNA solution, and colors correspond to different activity times  $t_a$  according to the color scale in (E). The dashed black curve is the control case in which no enzyme is added. After an initial increase compared to the control, the maximum plateau force  $F_p$  monotonically decreases with increasing activity time. (B) The time-course of  $F_p$  for solutions with 0 (open circles), 0.3 (filled cyan circles) and 0.5 (filled magenta circles)  $\text{U } \mu\text{g}^{-1}$  HaeIII. Each  $F_p$  value plotted is the average force value for  $x = 10\text{--}30 \text{ } \mu\text{m}$  with error bars representing standard error. The time-average of the  $0 \text{ U } \mu\text{g}^{-1}$  data shown is plotted as a dashed line for clarity. The initial increase in  $F_p$ , as well as the subsequent decrease, is slightly prolonged at lower enzyme concentration (cyan) but the functional form is preserved. (C) Force  $F$  shown in (A), normalized by the corresponding plateau force value  $F_p$  (light grey line) and plotted on a log-log scale shows that the distance at which the system yields to the terminal viscous-dominated plateau regime, *i.e.*, the yield distance  $x_y$ , decreases with increasing time  $t_a$ . The yield distance is quantified as the displacement at which  $F(x) = F_p e^{-1}$  (grey line). (D) The dependence of the yield distance on  $t_a$  shows similar trends as the plateau force shown in B. (E) Force relaxation  $F_R(t)$  measured during the relaxation phase, when the probe is held fixed at  $x = 30 \text{ } \mu\text{m}$ , normalized by the value at the beginning of the relaxation phase  $F_R(0)$ . A characteristic relaxation timescale  $\tau_0$  is quantified as the time at which  $F_R(t) = F_R(0)e^{-1}$  (grey line). (F) The trajectories of  $\tau_0$  versus  $t_a$  show similar trends as the plateau force and yield distance shown in B and D.

the lengthscale over which this elastic-like behavior crosses over to a strain-independent viscous plateau likewise decreases. We quantify this effective softening with increasing  $t_a$  by computing the yield distance  $x_y$  at which the force reaches  $F(x) = F_p/e$ , signifying yielding from elastic-like to viscous-dominated response. Rescaling each force trace by its corresponding plateau value  $F_p$ , clearly shows the decrease in yield distance, denoted by the  $x$  value at which each  $F(x)/F_p$  trace crosses the grey horizontal line in Fig. 2C, over the course of digestion. Moreover, the trajectory of  $x_y$  as a function of activity time (Fig. 2D) is highly similar to that of the plateau force (Fig. 2B), with an initial rise in  $x_y$ , followed by a sharp drop, before reaching a terminal minimum value at  $\sim 100$  min.

We expect this yield distance to be related to a characteristic confinement lengthscale that sets the elastic contribution to the response. For entangled linear polymers, this lengthscale is the diameter  $d_T$  of the confinement tube, formed by the surrounding polymers, that constrains polymer motion, resulting in elastic-like entropic stretching, rather than dissipative flow, in response to strain.<sup>23</sup> However, once an entangled linear

polymer has been strained a distance several times beyond  $d_T$  it can no longer maintain its entanglements with the surrounding polymers and instead the strain induces dissipative polymer disentanglement and flow.<sup>22</sup>

As described in the Introduction, because the tube model developed for linear polymers relies on the chains having free ends, the concept of ring entanglements is still debated and several models have been proposed to describe the dynamics of highly overlapping rings, in analogy to the reptation model. One such framework, the pom-pom ring (PPR) model, predicts an effective tube diameter for ring polymers that is smaller than its linear counterpart *via* the relation  $d_{T,R} \approx 0.7d_{T,L}$ ,<sup>44</sup> in line with the previously reported relation  $d_{T,R} \approx 0.75d_{T,L}$  for 115 kbp DNA.<sup>61</sup> These findings suggest that  $x_y$  should be proportionally larger when the solutions are primarily linear *versus* circular. Importantly, this description is only valid for systems with ample entanglements or other steric interactions, whereas unentangled polymers exhibit a very limited elastic regime (small  $x_y$ ) as the chains can more readily flow in response to strain.

The trajectory of  $x_y$  is qualitatively in line with this physical picture, as well as our interpretation of the plateau force data, as we expect initial linearization of the circular DNA to increase  $x_y$ , while fragmentation, which breaks and removes entanglements, to decrease  $x_y$  (Fig. 2D). Moreover, the maximum yield distance we measure, which we interpret to be when the polymers are mostly linear, is  $x_{y,L} \approx 480$  nm, compared to  $x_{y,0} \approx 315$  nm for the initial solution of circular polymers. The ratio  $x_{y,0}/x_{y,L} \approx 0.66$  is quite close to the theoretical and experimental ratios  $d_{T,R}/d_{T,L}$  discussed above, corroborating the concept that the yield distance is set by the tube diameter. Moreover, we can compute a tube diameter for linear DNA *via*  $d_{T,L} \approx (l_k l_e)^{1/2}$  where  $l_k \approx 100$  nm is the Kuhn length and  $l_e \approx L/Z \approx 200$  nm is the entanglement length,<sup>22,58,59</sup> providing  $d_{T,L} \approx 140$  nm and  $x_{y,L} \approx 3d_{T,L}$ . This relation quantitatively shows that yielding occurs once the polymers are strained beyond  $x \approx 3d_T$ , a reasonable assumption based on the physical picture described above. It follows that the ring polymers exhibit the same yielding behavior, with  $d_{T,0} \approx 98$  nm and  $x \approx 3d_T$ . Another important constraint lengthscale that likely contributes to the initial increase in elastic-like response is that associated with the threaded rings that we expect to be present at early digestion times. Threaded rings are confined laterally by a distance  $\sim R_{G,R} \approx 135$  nm which is quite close to  $d_{T,L} \approx 140$  nm. So, while we cannot parse the contributions of tube confinement *versus* lateral pinning from the yielding behavior, the observed stress overshoots suggest that both mechanisms play a role. Finally, the large drop in  $x_y$  as digestion proceeds is consistent with fragmentation that reduces and eventually eliminates entanglements.

Next, we turn to the relaxation of the force following straining to determine how the longest relaxation timescales vary during digestion. As shown in Fig. 2E, the force relaxation follows similar trends as the response during strain, with the decay in  $F_R$  first becoming slower, seen as a shallower slope for  $t_a = 10$  min, followed by rapidly increasing decay rates, seen as steeper slopes that approach a relative constant value for  $t_a > 50$  min. To approximate a single characteristic relaxation time  $\tau_0$  for the system at each point during digestion, we compute the measurement time at which  $F_R(t) = F_R(0)/e$  (Fig. 2E). As shown in Fig. 2F,  $\tau_0$  exhibits the same dependence on activity time as  $F_p$  and  $x_y$ , corroborating the physical picture that the elastic contribution to the nonlinear stress response initially increases, due to linearization of the DNA circles that increases entanglements and threadings, after which yielding to a viscous-dominated regime ensues due to fragmentation that reduces steric constraints.


To shed light on the principal mechanism underlying the force relaxation, we note that previous studies examining the nonlinear response of entangled DNA solutions have shown that the relaxation modes that generally contribute to the force relaxation are elastic Rouse relaxation, which occurs over the Rouse time  $\tau_R$ ; reptation, which occurs over the disengagement time  $\tau_D \approx 3Z\tau_R$ ; and entanglement, which occurs over the entanglement time  $\tau_e \approx Z^{-2}\tau_R$ .<sup>22,48,57,61,66,67</sup> As discussed above, the Rouse times for linear, ring and supercoiled topologies are

$\tau_{R,L} \approx 136$  ms,  $\tau_{R,R} \approx 42$  ms, and  $\tau_{R,S} \approx 19$  ms. For linear DNA, the disengagement time is  $\tau_{D,L} \approx 3245$  ms, and for rings we use the relation  $\tau_{D,R} \approx 0.3\tau_{D,L}$ , predicted by the PPR model and experimentally validated for DNA, to estimate  $\tau_{D,R} \approx 973$  ms.<sup>44,48,57,61</sup> While there is no analogous relation for  $\tau_{D,S}$  for supercoiled constructs, nor is it clear the extent to which they participate in entanglements, we can roughly estimate  $\tau_{D,S}/\tau_{D,L} \approx \tau_{R,S}/\tau_{R,L}$ , which holds for ring polymers, to approximate  $\tau_{D,S} \approx 450$  ms. To provide a rough estimate of topology-dependent  $\tau_e$  values, we compute  $\tau_{e,L} \approx 1.4$  ms, and use the scaling relation  $d_T^{-2} \sim Z$  and  $d_{T,R} \approx 0.7d_{T,L}$  to approximate  $Z_R \approx 0.5Z_L$ , from which it follows that  $\tau_{e,R} \approx 1.6$  ms. If entanglements play a role in the relaxation of supercoiled constructs, then it is reasonable to assume  $\tau_{e,S}$  is of similar order of magnitude as  $\tau_{e,L}$  and  $\tau_{e,R}$ , given the similarity between  $\tau_{e,R}$  and  $\tau_{e,L}$ .

Considering these values, it is clear that reptation is not a primary mode of nonlinear stress relaxation, as all  $\tau_D$  values are  $\sim 1$ – $2$  orders of magnitude longer than  $\tau_0 \approx 6$ – $30$  ms (Fig. 2F). Our estimated  $\tau_e$  values are likewise shorter than  $\tau_0$ , suggesting minimal contribution. Our measured  $\tau_0$  values are, however, comparable to  $\tau_R$  for the different topologies, suggesting that relaxation proceeds primarily *via* Rouse modes. This interpretation is consistent with our finding that all solutions reach a terminal viscous plateau in which entanglements do not play an appreciable role, and with previous studies that have shown that in certain cases even well-entangled linear DNA subject to nonlinear strains can completely relax their stress *via* Rouse relaxation.<sup>64</sup> Given the topological complexity and polydispersity of our solutions during digestion, we do not expect exact agreement with predictions for monodisperse systems, but order of magnitude agreement with predictions serves as a guide towards underlying mechanisms.

### Linear microrheological response shows opposite dependence on activity time at low frequencies *versus* high frequencies

To determine how topological digestion of the DNA solutions alters the linear viscoelastic response, we measure the thermal fluctuations of a trapped probe to extract the frequency-dependent linear elastic and viscous moduli  $G'(\omega)$  and  $G''(\omega)$  (Fig. S3, ESI<sup>†</sup>). These measurements probe smaller lengthscales and slower strain rates than nonlinear measurements, with a range of accessible frequencies of  $\omega = 1$ – $300$  rad  $s^{-1}$ , corresponding to timescales of  $t = 2\pi\omega^{-1} = 21$  ms– $6.25$  s, which are all longer than the nominal nonlinear strain timescale  $t \approx \dot{\gamma}^{-1} \approx 17$  ms. To examine the elastic contribution to the linear stress response, we evaluate  $g'(\omega) = G'(\omega)/|G^*(\omega)|$  where  $|G^*(\omega)| = [(G'(\omega))^2 + (G''(\omega))^2]^{1/2}$  is the magnitude of the complex modulus, and  $g'(\omega)$  ranges between 1 for a purely elastic response to 0 for purely viscous response. Fig. 3A reveals distinct dependence on the activity time at high frequencies *versus* low frequencies. Specifically, for  $\omega = 30$ – $300$  rad  $s^{-1}$ ,  $g'(\omega)$  exhibits a similar dependence on  $t_a$  as the nonlinear response metrics, with  $g'(\omega)$  generally decreasing with increasing  $t_a$ . However, for much lower frequencies of  $\omega = 1$ – $5$  rad  $s^{-1}$ , we observe a counterintuitive nearly opposite trend, with  $g'(\omega)$  first



**Fig. 3** Linear OTM reveals distinct trajectories of viscoelastic moduli at high and low frequencies during topological digestion. (A) Fractional elasticity  $g'(\omega) = \frac{G'(\omega)}{|G^*(\omega)|}$  versus frequency  $\omega$  for increasing activity times  $t_a$  after adding  $0.5 \text{ U } \mu\text{g}^{-1}$  of Haelll, denoted by the color scale. Control data ( $0 \text{ U } \mu\text{g}^{-1}$ ) is shown as open symbols. The dependence of  $g'(\omega)$  on activity time  $t_a$  differs at low frequencies (blue,  $\omega = 1\text{--}5 \text{ rad s}^{-1}$ ) compared to the high frequency regime (yellow,  $\omega = 30\text{--}300 \text{ rad s}^{-1}$ ). (B) Trajectories of  $g'$ , averaged over all frequencies in the high frequency regime  $[(g')(t_a)]_H$  (top, yellow) and low frequency regime  $[(g')(t_a)]_L$  (bottom, blue) for Haelll:DNA stoichiometries of  $0.3 \text{ U } \mu\text{g}^{-1}$  (filled cyan circles) and  $0.5 \text{ U } \mu\text{g}^{-1}$  (filled magenta circles)  $\text{U } \mu\text{g}^{-1}$ , compared to the average value for the undigested control case (dashed line). (C) Complex viscosity  $\eta^*(\omega)$  versus frequency  $\omega$  for increasing activity times  $t_a$  after adding  $0.5 \text{ U } \mu\text{g}^{-1}$  of Haelll, denoted by the color scale. Control data ( $0 \text{ U } \mu\text{g}^{-1}$ ) is shown as open symbols. The power-law scaling  $\eta^*(\omega) \sim \omega^{-\beta}$  is generally steeper (larger  $\beta$ ) at low frequencies (blue region) compared to high frequencies (yellow region), and depends on  $t_a$ . (D) Trajectories of  $\beta(t_a)_H$  (top, yellow) and  $\beta(t_a)_L$  (bottom, blue) determined via fits of  $\eta^*(\omega)$  to the high and low frequency ranges, respectively, at each activity time  $t_a$  for Haelll:DNA stoichiometries of  $0.3 \text{ U } \mu\text{g}^{-1}$  (filled cyan circles) and  $0.5 \text{ U } \mu\text{g}^{-1}$  (filled magenta circles)  $\text{U } \mu\text{g}^{-1}$ , compared to the average value for the undigested control case (dashed line). At high frequencies the time dependence matches that observed in the nonlinear regime, while a nearly opposite trend is seen at low frequencies.

decreasing below that of the control, followed by an increase, suggestive of increased elasticity upon fragmentation.

By averaging  $g'(\omega)$  over all frequencies in each regime, we can more clearly compare trajectories of  $g'$  in the high and low frequency ranges,  $[(g')(t_a)]_H$  and  $[(g')(t_a)]_L$  (Fig. 3B). The difference from the control data is less extreme compared to the nonlinear regime metrics (Fig. 2), but the different trends are obvious.  $[(g')(t_a)]_H$  follows a similar trajectory as the nonlinear regime data, reaching a  $t_a$ -independent plateau that is lower than the control. Conversely, the terminal  $t_a$ -independent plateau in  $[(g')(t_a)]_L$  is higher than the control. The takeaway is that, at long times, the solutions appear to be locally more sterically constrained and elastic-like following digestion, whereas, at short times, the local response is suggestive of fewer constraints and more dissipative contributions.

We also evaluate the complex viscosity  $\eta^*(\omega) = |G^*(\omega)|\omega^{-1}$  which shows shear-thinning  $\eta^*(\omega) \sim \omega^{-\beta}$  for all frequencies and activity times. We previously reported scaling exponents of  $\beta \approx 0.6$  and  $\beta \approx 0.4$  for concentrated solutions of linear and ring DNA, as well as a higher value,  $\beta \approx 0.7$  for ring-linear

blends,<sup>48</sup> whereas we found  $\beta \approx 0.5$  for unentangled semidilute polymers and  $\beta \rightarrow 0$  as chains no longer overlap ( $c < c^*$ ).<sup>68</sup> Fig. 3c reveals that the thinning exponents are generally higher in the low frequency range ( $\omega = 1\text{--}5 \text{ rad s}^{-1}$ ) versus the high frequency regime ( $\omega = 30\text{--}300 \text{ rad s}^{-1}$ ), indicative of increased constraints (*i.e.*, entanglements, threadings). More importantly, the trajectory of  $\beta$  with activity time displays the same switch between high and low frequency regimes. Namely, at low frequencies,  $\beta(t_a)_L$  generally increases with activity time, similar to  $[(g')(t_a)]_L$ , reaching a plateau value that is higher than that of the control. Conversely,  $\beta(t_a)_H$  generally decreases with increasing  $t_a$ , with the trajectory mirroring  $[(g')(t_a)]_H$ , reaching a plateau that has a lower value than the control. These results indicate that the local dynamics are dominated by different mechanisms at short times  $t = 2\pi\omega^{-1}$  of  $21 \text{ ms} \lesssim t_s \lesssim 210 \text{ ms}$  versus long times  $1.25 \text{ s} \lesssim t_l \lesssim 6.3 \text{ s}$ .

To consider the differences between these timescales relative to the intrinsic relaxation timescales of the system, we recall the predicted topology-dependent values for  $\tau_D$ ,  $\tau_R$  and  $\tau_e$  discussed in the previous section. The long-time regime

extends to times longer than all of the predicted relaxation timescales, with  $\tau_{D,L} \simeq 3.2$  s falling within the long-time range and  $\tau_{D,R} \simeq 0.97$  s lying just short of  $t_s$ , so we expect reptation of the long chains to dominate the response. Conversely, the short-time regime is below  $\tau_D$  and above  $\tau_e$  for all topologies, suggesting that unentangled Rouse dynamics likely dictate the stress response and that the primary contributors to this short-time stress response are ostensibly unentangled circular DNA and shorter linear fragments that have terminal relaxation times much faster than  $\tau_D$ . We will consider this argument further in the Discussion.

### DNA transport exhibits counterintuitive slowing during enzymatic fragmentation

To shed light on the timescale-dependent impact of topological alteration on rheology, we turn to measuring how digestion alters the apparent diffusion of DNA molecules in the topologically-active solutions at long timescales. We use DDM for this purpose, rather than single-molecule tracking, to allow for more statistics in the long-time regime (Fig. 1G–J and 4). As shown in Fig. 4A, DDM is able to extract statistically meaningful intermediate scattering functions (ISFs) for lag times up to  $t = 10$  s. Visual inspection of the ISFs for different activity times  $t_a$  (Fig. 4A) reveals that ISFs decay more slowly as  $t_a$  increases, indicative of slowing down of dynamics during digestion, in qualitative alignment with the low-frequency microrheology data (Fig. 3). By fitting the ISFs for each wave vector  $q = 2\pi\lambda^{-1}$  to an exponential (see Methods) we extract the characteristic decay time  $\tau(q)$  for each time  $t_a$ , which is a measure of how quickly density fluctuations at a lengthscale  $\lambda$  become decorrelated.<sup>69–71</sup> While we expect heterogeneous ensembles of molecules at each stage of the digestion process, which will have different intrinsic relaxation timescales, the separation between the various topology-dependent and size-dependent timescales, as well as the likely broad spectrum of values, limits the ability to resolve these distinct timescales using DDM. Rather, our measurements provide an average decay time that

reflects the collective fluctuations of the ensemble of molecules. Increasing  $\tau(q)$  values at a given  $q$  indicate slowing of dynamics, which Fig. 4B shows occurs as  $t_a$  increases. Moreover all  $\tau(q)$  data roughly follow power-law scaling  $\tau(q) \sim q^{-2}$ , expected for diffusive dynamics (rather than subdiffusive or superdiffusive), and the diffusion coefficient is computed *via*  $D = \tau^{-1}q^{-2}$ .<sup>35,70,72,73</sup>

Fig. 4C plots  $D$  versus  $t_a$ , as determined from fits of  $\tau(q)$ , revealing a strikingly similar trajectory as  $[(g')(t_a)]_L$  and  $\beta(t_a)_L$ , with an initial slowing of dynamics (expected for increased elasticity and steric constraints) until reaching a plateau at  $t_a \approx 100$  s for  $0.5 \text{ U } \mu\text{g}^{-1}$ . The lower stoichiometry ( $0.3 \text{ U } \mu\text{g}^{-1}$ ) decays more slowly due to slower digestion, as seen in the OTM data, albeit to a lesser extent.

To determine how general this result is, we perform DDM on solutions undergoing digestion by another enzyme, ApoI, which cuts the DNA at a different number of recognition sites compared to HaeIII (8 *vs.* 17). The resulting fragments also have sticky overhangs rather than the blunt ends that HaeIII cleaving provides (Fig. S2, ESI†). As shown in Fig. 5A,  $\tau(q)$  increases with increasing activity time, just as we see for HaeIII (Fig. 4B), and follows a trajectory similar to that of  $0.5 \text{ U } \mu\text{g}^{-1}$  HaeIII. To ensure that this slowing is not an artefact of our DDM analysis, we also track single DNA molecules in the topologically-active solutions and evaluate their mean-squared displacements (MSD) *versus* lag time. Fig. 5B shows that indeed the diffusive slowing is preserved in the tracking data, with the MSDs decreasing with increasing activity time, indicative of slower dynamics, and the motion is primarily diffusive, as indicated by the linear scaling of the MSDs with lag time (*i.e.*,  $\text{MSD} \sim t^\alpha$  with  $\alpha \simeq 1$ ). From linear fits to the MSDs, we compute diffusion coefficients, which correlate remarkably well with those measured *via* DDM over the full time-course of activity (Fig. 5C). The takeaway: for different measurement techniques, enzymes and stoichiometries, the slowing of DNA dynamics is universal and trajectories are largely insensitive to the details.





Fig. 4 DDM of fluorescent-labeled DNA shows that HaeIII-mediated fragmentation of DNA surprisingly slows the diffusion of DNA. (A) Intermediate scattering function (ISF) of topologically-active DNA solution as a function of lag time  $\Delta t$  at different activity times  $t_a$  during digestion by  $0.5 \text{ U } \mu\text{g}^{-1}$  HaeIII, color coded according to the legend. As digestion proceeds (as  $t_a$  increases) the time at which ISFs decay is extended, indicating slowing of dynamics. (B) Decay time  $\tau(q)$  vs. wave vector  $q$ , computed *via* fitting the ISF, shows diffusive dynamics for all times  $t_a$ , *i.e.*,  $\tau(q) \sim q^{-2}$  but the value of  $\tau(q)$  at each  $q$  increases with increasing  $t_a$ , indicating lower diffusion coefficients  $D$  (slower dynamics), which can be computed *via* the relation  $\tau(q) = 1/Dq^2$ . (C) Diffusion coefficients  $D$ , determined from  $\tau(q)$  *versus* activity time for HaeIII:DNA stoichiometries of  $0.3 \text{ U } \mu\text{g}^{-1}$  (filled cyan circles),  $0.5 \text{ U } \mu\text{g}^{-1}$  (filled magenta circles) and 0 (open circles, dashed line denotes average).





**Fig. 5** The apparent slowing of DNA transport observed at long times is universal across analysis methods as well as enzyme type and concentration. (A) Decay time  $\tau(q)$  vs. wave vector,  $q$  computed identically as in Fig. 4B, for DNA solutions undergoing digestion by the multi-site restriction endonuclease, Apol, at a Apol:DNA stoichiometry of  $0.3 \text{ U } \mu\text{g}^{-1}$ . Apol digestion produces 8 fragments with sticky overhangs, compared to 17 blunt-ended fragments produced by HaellI. The observed  $q$ -scaling of  $\tau(q)$  and its increase with increasing  $t_a$  are very similar to those seen in Fig. 4B for  $0.5 \text{ U } \mu\text{g}^{-1}$  HaellI. (B) Tracking of single fluorescent-labeled DNA molecules in the same conditions shown in Fig. 4A and B ( $0.5 \text{ U } \mu\text{g}^{-1}$  HaellI) shows that the mean-squared displacement (MSD) of the molecules are largely diffusive, *i.e.*,  $\text{MSD} \sim t^\alpha$  with  $\alpha = 1$ , similar to results from DDM. Moreover, the magnitude of the MSD at any given lag time  $t$  decreases with increasing activity time  $t_a$ , indicating a slowing of molecular diffusion, similar to that seen in A and Fig. 4. (C) Diffusion coefficients  $D$  calculated using DDM and single-molecule tracking (SMT) for varying concentrations ( $0.3 \text{ U } \mu\text{g}^{-1}$  (cyan, teal),  $0.3 \text{ U } \mu\text{g}^{-1}$  (magenta, pink)) and types (HaellI (H), Apol) of enzymes show universal slowing of diffusivity during digestion that is rapid initially and then weaker, appearing to asymptote to a steady-state value when digestion is complete. Apol and SMT data are determined from A and B while HaellI DDM data are the same as in Fig. 4C.


## Discussion

The data discussed above shows stark differences between how linearization and fragmentation alter the rheology and dynamics at different timescales. The short-time dynamics make intuitive sense. Namely, as the DNA gets chopped into smaller and smaller pieces, we expect the elastic contribution to the stress response to go down. However, upon initial linearization of the rings this elastic contribution increases as the degree of polymer overlap and entanglement density increase and threading events become possible. These general trends were also observed in recent particle-tracking micro-rheology experiments that showed that linearization of circular DNA lead to a steady increase in viscosity, while fragmentation of linear DNA led to a steady decrease.<sup>13</sup>

To validate this physical picture and quantitatively map the topological composition of the solution to the rheological properties during digestion, we perform time-resolved gel electrophoresis experiments to determine the fraction of ring, supercoiled, full-length linear, and fragmented linear DNA at each measurement time point during digestion (Fig. 6A and B). The key features that align with our interpretation of the rheology are: (1) linearization begins immediately and the fraction of linear chains steadily increases over the first  $\sim 50$  min before reaching a maximum value of  $\sim 50\%$  that persists until  $\sim 125$  min after which the fraction of full-length linear DNA decreases back to  $0\%$ ; (2) rings persist until  $\sim 100$  min, even as fragmentation begins; (3) supercoiled DNA is digested more rapidly than rings and are completely eliminated within  $\sim 50$  min. These observations suggest that the initial sharp increase in elasticity is due to linearization that not only increases the entanglement density but also increases the number of threading events due to the persistent presence of rings. As the fraction of linear chains increases further and the fraction of rings likewise decreases, the number of threading

events decreases, which causes the elastic-like contributions to decrease despite the increased entanglement density. The subsequent slower decrease in elasticity from  $\sim 50$  min to  $\sim 125$  min aligns with the onset of fragmentation and elimination of circular DNA while maintaining  $\sim 50\%$  full-length DNA. Here, circular DNA molecules are, on average, exchanged for linear fragments, which decreases threading events and entanglements. After  $\sim 125$  min, the fraction of fragments exceeds that of full-length chains, such that there are insufficient entanglements to support elastic-like stress response. Once this broken percolation of entanglements is reached, further fragmentation has little impact on the rheology. We also note that the rapid digestion of supercoiled molecules supports our conjecture that supercoiled constructs do not play a significant role in the relaxation dynamics of the digesting solution.

To understand the long-time behavior, we turn to recent work, described in the Introduction, that examined the bulk rheology of topologically-active solutions of 11 kbp circular DNA and their composites with dextran.<sup>14</sup> In these studies, measurements were performed at  $\omega = 1 \text{ rad s}^{-1}$ , which is the lowest frequency (longest time) that we evaluate here, over the course of digestion by a multi-site cutting enzyme that cleaved the rings into 13 fragments of different lengths. These experiments showed that fragmentation actually led to an abrupt increase in  $G'$  which only occurred following complete digestion (rather than initial linearization). Similar abrupt jumps in  $G'$  were observed in composites with dextran, but the onset was much earlier in the digestion process. This early onset was shown to be due to crowding of the DNA by dextran that facilitated clustering of the DNA, *via* depletion interactions, that increased the local polymer overlap and entanglements within the clusters. These 'slow' clusters, where the DNA is more constrained and diffuses more slowly, dominate the system dynamics at long times, but are frozen



**Fig. 6** Different relaxation mechanisms manifest at short and long timescales during digestion of circular DNA into polydisperse fragments. (A) Time-resolved gel electrophoresis of DNA undergoing HaeIII digestion. Each lane shows the topological state of the DNA at a given time  $t_a$ , in 10 min intervals, as listed in min on the x-axis. Far-left lane is the  $\lambda$ -HindIII molecular weight marker with the size of linear DNA corresponding to select bands listed in basepairs (bp). Translucent boxes with increasingly lighter shades of green denote supercoiled, ring, full-length linear and fragmented linear DNA, respectively. Vertical dashed lines indicate  $t_a = 50, 100, 150, 200$  min. (B) Mass fraction of supercoiled (dark green triangles), ring (olive circles), full-length linear (green squares) and fragmented linear (light green x's) DNA at each time  $t_a$  during digestion, determined *via* band intensity analysis of the gel shown in A. (C) Various metrics of confinement  $X$ , shown in previous figures, normalized by their  $t_a = 0$  values  $X_{t_a=0}$ , and averaged over the 0.3 and 0.5 U  $\mu\text{g}^{-1}$  HaeIII conditions, are color coded according to the timescale that the measurement probes (short/yellow or long/blue). The different symbols and shades of each color indicate the metric or measurement technique according to the legend. For DDM and nonlinear OTM, the normalized effective viscosities  $\eta/\eta_{t_a=0}$  are plotted, computed as  $\eta/\eta_{t_a=0} = D_{t_a=0}/D$  and  $\eta/\eta_{t_a=0} = F_p/F_{p,t_a=0}$ , respectively. Short timescale metrics universally indicate that topological activity increases dissipation and mobility, while the opposite is observed at long timescales. (D) Three decades of timescales that different measurement techniques and metrics probe, colorized according to short (yellow shades) and long (blue shades) timescales. Low-frequency linear OTM, DDM, and SMT probe long timescales ( $\sim 0.2$ – $6$  s) while high-frequency linear OTM and nonlinear OTM probe fast timescales ( $\sim 17$ – $200$  ms).

out at short times, below the relaxation timescales of the clustered DNA.

Similar effects of DNA crowding by dextran have been seen in bulk rheology measurements of steady-state DNA-dextran composites (no digestion).<sup>74</sup> In the absence of dextran, clustering of longer DNA chains is still expected to occur *via* crowding by the smaller fragments that have similar radius of gyration as dextran.<sup>14</sup> Specifically, enhanced self-association and clustering of molecules in the presence of crowders occurs as the crowding population seeks to maximize its entropy by maximizing its available volume, driving larger molecules together. This depletion effect only occurs, *i.e.*, is only entropically favorable, when the crowders are sufficiently small and at high enough concentration, as compared to molecules being crowded. For reference, the radius of gyration of the dextran used in the previous studies was  $R_G \approx 19$  nm compared to  $R_G \approx 280$  nm and  $R_G \approx 196$  nm for the two DNA lengths studied. Here,  $R_G \approx 103$ – $213$  nm for DNA (depending on the topology), while the radius of gyration for the average fragment length at complete digestion,  $\langle l_f \rangle \approx 115$  nm, is  $R_G \approx 20$  nm and for the smallest length fragment  $l_{f,\text{min}} \approx 6.3$  nm (19 basepairs) is  $R_G \approx 5$  nm. Given the similar sizes of dextran and the DNA fragments in our current work and previous studies, it is reasonable to assume that similar depletion interactions may be at play. Moreover, the roughly order of magnitude difference between the size of DNA and its crowders (fragments or dextran) in all cases, suggests that we may expect similar mechanisms to drive the non-equilibrium rheological behavior.

Based on this knowledge, we hypothesize that the increased elasticity and slowed diffusion that fragmentation elicits at long timescales is a result of the depletion-driven formation of local clusters of long DNA *via* crowding by shorter fragments. The DNA molecules within these clusters are more entangled and slower than those in the surrounding solution, and it is the dynamics of these slow clusters that dominate the rheology and diffusion of the system at long times, *i.e.*, those comparable to the relaxation of the clustered DNA. This slowing is not seen at short times as the relaxation of these slow molecules is frozen out on this fast measurement timescale. To corroborate this physical picture, we point to the larger relative drop in DNA diffusion coefficients, measured *via* DDM and single-molecule tracking, compared to the relative increase in  $\langle g' \rangle$  and  $\langle \beta \rangle$  measured *via* OTM (Fig. 6C). The former measurements rely on fluorescent-labeled DNA to report dynamics while the latter does not. However, DNA below  $\sim 2$  kbp cannot be seen in typical fluorescence microscopes due to resolution limits and the inability to add enough fluorophores to smaller chains. Nevertheless, we were still able to visualize labeled DNA throughout the duration of digestion (Fig. 1G), suggesting that the DNA that we were tracking were the longest fragments, while the shorter ones were unable to contribute to the signal. Thus, we expect our DDM and tracking data to be representative of only the slow-moving clusters of longer fragments. OTM measurements on the other hand measure the local rheology in the vicinity of the trapped bead, which could be in a cluster or in the faster surroundings, so the apparent increase in elasticity is more subdued.

To more directly compare the changes in rheological properties during digestion, measured over different timescales and with different techniques (Fig. 6D), we plot together the trajectories of all metrics shown in Fig. 2–4, normalized by their  $t_a = 0$  value, and averaged over the two stoichiometries (Fig. 6C). For more straightforward comparison of metrics that increase with increased elasticity and slowing of dynamics, we convert our DDM-measured normalized diffusion coefficients  $D(t_a)/D_{t_a=0}$  to viscosity *via* the relation  $D \sim \eta^{-1}$ , and our maximum plateau force  $F_p/F_{p,t_a=0}$  to viscosity *via* the relation  $F_p = \eta v$ . Fig. 6C shows the trajectories of these long-time (DDM) and short-time (nonlinear OTM) viscosities, along with  $\langle g' \rangle$  and  $\langle \beta \rangle$  for both short and long measurement timescales. The short-time trajectories of  $\eta(t_a)/\eta_{t_a=0}$ ,  $\langle g' \rangle/\langle g' \rangle_{t_a=0}$ , and  $\langle \beta \rangle/\langle \beta \rangle_{t_a=0}$ , are similar for all three metrics, with the nonlinear OTM viscosity decaying to a slightly lower value, presumably because of the faster timescales probed ( $\sim 18$  ms *vs.* 21–210 ms). The long-time trajectories show a more pronounced difference, with  $\eta(t_a)/\eta_{t_a=0}$  increasing by  $\sim 80\%$  over the course of digestion, while  $\langle g' \rangle/\langle g' \rangle_{t_a=0}$ , and  $\langle \beta \rangle/\langle \beta \rangle_{t_a=0}$  increase by  $\sim 22\%$ . This significant difference corroborates the effect described above of only being able to visualize the larger fragments in DDM (as well as single-molecule tracking, Fig. 5B).

## Conclusion

In conclusion, we have presented the results of comprehensive measurements of concentrated circular DNA pushed out-of-equilibrium *via* the action of multi-site cleaving enzymes that linearize and chop up DNA circles. Using linear and nonlinear optical tweezers microrheology, as well as DDM and SMT, we span over 3 orders of magnitude in measurement timescales to reveal a surprising timescale-dependent variation of viscoelastic moduli and DNA diffusivity during enzymatic activity. Namely, at short timescales, well below the nominal disengagement times for full-length circular and linear DNA, probed *via* nonlinear OTM (Fig. 2) and high-frequency linear OTM (Fig. 3), digestion leads to an initial increase in metrics of elasticity and confinement, *i.e.*,  $g'$ ,  $\beta$ ,  $\eta$ ,  $x_p$ ,  $\tau_0$ , due to linearization of the circular DNA that increases entanglements and threading events. This increase is followed by a rapid decrease as fragmentation continues, due to disruption of entanglements and decreased polymer overlap, after which a plateau is reached once fragments dominate the solution composition. At long timescales, comparable to or longer than the disengagement times for the full-length DNA, probed *via* low-frequency OTM (Fig. 3), DDM (Fig. 4) and SMT (Fig. 5), we observe remarkably opposite trajectories of metrics, whereby fragmentation causes a rapid increase in elasticity metrics, *i.e.*,  $g'$ ,  $\eta$ ,  $\beta$ , after an initial transient drop, that stabilizes to a plateau value at nearly the same activity time as the short-time metrics.

We conjecture that these results arise from the different relaxation mechanisms and molecular lengthscales that the two regimes can detect, along with crowding induced self-association of longer DNA chains by short fragments. At short

timescales, the relaxation modes of the long DNA strands are prohibitively slow to be detected and are thus frozen out, leaving the shorter fragments and more weakly entangled molecules to dominate the signal. As fragmentation proceeds, these shorter chains become shorter, which decreases overlap and entanglements that provide elasticity. Moreover, the depletion of the longer chains by the shorter chains increases the volume available to the shorter chains, making them effectively more dilute and mobile, further decreasing any elastic contributions to the stress response. Conversely, at long lengthscales, the motion of the slow-moving clusters of long DNA dominates the response. As fragmentation proceeds, the depletion interaction strength increases as more and more fragments are created, increasing the self-association of longer chains. This effect increases the local entanglement density and constraints, leading to an increase in elasticity metrics and slowing of diffusion. While this physical picture is, at present, a hypothesis for which we do not have unequivocal evidence, we hope that our intriguing results will spur new theoretical and simulation studies that can validate and shed further light on the proposed mechanisms.

Our results highlight the ability to use topology to drive non-equilibrium restructuring and viscoelastic alterations that can be tuned at the molecular level to elicit starkly different trajectories at different spatiotemporal scales. These hierarchical reconfigurable materials may prove useful in applications in sequestration, self-healing, filtration, and micro-actuation. More generally, the use of enzymes to impart activity into steady-state systems offers a mechanism for self-assembly, reconfigurability, and self-healing with timing and strength that can be precisely programmed by the stoichiometries and types of enzymes and substrate.

## Methods

### DNA

We prepare solutions of double-stranded DNA, of length  $L = 5.9$  kilobasepairs (kbp), *via* replication of pYES2 plasmid constructs in *Escherichia coli* followed by extraction, purification and concentrating as described previously.<sup>75–77</sup> Briefly, cultures of *E. coli* cells containing pYES2 are grown from frozen glycerol stocks, after which the DNA is extracted from cells *via* alkaline cell lysis, renatured *via* treatment with an acidic detergent, precipitated *via* isopropanol, and resuspended in nanopure deionized water (DI). To remove RNA and proteins the solution is treated with Rnase A followed by phenol-chloroform extraction and dialysis. The purified DNA solution is further concentrated *via* rotary vacuum concentration and stored at 4 °C.

We use gel electrophoresis and band intensity analysis, employing Life Technologies E-Gel Imager and Gel Quant Express software, to determine a stock DNA concentration of  $c = 8$  mg mL<sup>-1</sup> and a topological distribution of  $\sim 62\%$  relaxed circular (ring) DNA and  $\sim 38\%$  supercoiled circular DNA. The radius of gyration for each topology is  $R_{G,R} \simeq 135$  nm and  $R_{G,S} \simeq 103$  nm, from which we compute an overlap

concentration of  $c^* \simeq \left[ \frac{3M}{4\pi N_A} \right] [0.55R_{G,R}^3 + 0.45R_{G,S}^3]^{-1} \simeq 750 \mu\text{g mL}^{-1}$ .<sup>21</sup> For reference, for a solution of 100% linear DNA at the same concentration, the quantities are  $R_{G,L} \simeq 213 \text{ nm}$  and  $c^* \simeq 157 \mu\text{g mL}^{-1}$ .

### Restriction endonucleases

We use the high-fidelity restriction endonuclease HaeIII (New England BioLabs) to introduce topological activity into the DNA solution. HaeIII cleaves pYES2 at 17 different 5'-GG/CC-3' recognition sites (Fig. 1a), converting each circular DNA to linear topology (the first cut) and then cutting the resulting linear strand into 17 blunt-ended linear fragments of different lengths ranging from  $l_{f,\text{min}} = 19 \text{ bp}$  to  $l_{f,\text{max}} = 1642 \text{ bp}$ , with an average fragment length of  $l_f \simeq 345 \text{ bp}$  (Fig. S1, ESI†). For all data shown we use a HaeIII:DNA stoichiometry of  $0.3 \text{ U } \mu\text{g}^{-1}$  or  $0.5 \text{ U } \mu\text{g}^{-1}$ , chosen to be low enough to ensure that any topological activity can be considered frozen out on the time-scale of each measurement ( $\sim 10 \text{ s}$ ) but fast enough that complete digestion is achieved over the course of each 4 hour experiment. We also performed a subset of experiments using the high-fidelity restriction endonuclease ApoI (New England BioLabs) at stoichiometry of  $0.3 \text{ U } \mu\text{g}^{-1}$  DNA. ApoI cleaves pYES2 at 8 different 5'-R/AATY-3' recognition sites, resulting in 8 fragments with sticky overhangs and lengths from  $l_{f,\text{min}} = 11 \text{ bp}$  to  $l_{f,\text{max}} = 2585 \text{ bp}$ , with  $l_f \simeq 738 \text{ bp}$  (Fig. S2, ESI†).

### Sample preparation

To prepare topologically active DNA solutions for all experiments, we dilute the stock DNA to a final concentration of  $6 \text{ mg mL}^{-1}$  with 1/10th volume of  $10\times$  CutSmart Buffer ( $0.5 \text{ M}$  Potassium Acetate,  $0.2 \text{ M}$  Tris-acetate,  $0.1 \text{ M}$  Magnesium Acetate; New England BioLabs),  $0.1\%$  Tween (to prevent non-specific interactions), and the appropriate volume of HaeIII. We mix the solution by pipetting up and down with a wide-bore pipet tip, and add the enzyme last, noting the time at which it is added, which we define as  $t_a = 0$ . The buffer conditions and temperature ( $20 \text{ }^\circ\text{C}$ ) provide good solvent conditions for the DNA.<sup>58,65,76,78,79</sup>

For DDM and single-molecule tracking, we include  $2\text{--}20 \mu\text{g mL}^{-1}$  of fluorescent-labelled DNA and an oxygen scavenging system ( $45 \mu\text{g mL}^{-1}$  glucose,  $43 \mu\text{g mL}^{-1}$  glucose oxidase,  $7 \mu\text{g mL}^{-1}$  catalase,  $5 \mu\text{g mL}^{-1}$  beta-mercaptoethanol) to inhibit photobleaching. The DNA is labelled with covalent dye Mirus-488 (Mirus Bio) at a dye : basepair ratio  $1 : 5$  using the Mirus Label IT Nucleic Acid Labeling Kit (Mirus Bio).

For optical tweezers microrheology measurements we add a trace amount of  $4.5 \mu\text{m}$  diameter polystyrene microspheres (Polysciences) that are coated with BSA to inhibit non-specific interactions with the DNA and enzymes.

We load the final  $10 \mu\text{L}$  sample into a microscope sample chamber created by fusing together a glass coverslip and microscope slide with a heated parafilm spacer, and seal the chamber with UV glue. The chamber surfaces are passivated

with  $10 \text{ mg mL}^{-1}$  BSA to prevent non-specific adsorption of the DNA or enzyme.

For time-resolved gel electrophoresis experiments, we increase sample volume to  $40 \mu\text{L}$  and keep the solution in a microcentrifuge tube during digestion.

### Optical tweezers microrheology

We use optical tweezers microrheology to determine the linear and nonlinear rheological properties over the course of digestion (Fig. 1B–F). The optical trap, built around an Olympus IX71 epifluorescence microscope, is formed from a  $1064 \text{ nm}$  Nd:YAG fiber laser (Manlight) focused with a  $60\times$   $1.4\text{NA}$  objective (Olympus), as described previously.<sup>63,64,68</sup> The force exerted by the DNA on a trapped bead is determined by recording the laser beam deflections *via* a position sensing detector (Pacific Silicon Sensors) at  $20 \text{ kHz}$ . The trap is calibrated for force measurement using Stokes drag method.<sup>63,64,68</sup> All microrheological data is recorded at  $20 \text{ kHz}$ , and each measurement is taken with a new microsphere in an unperturbed location. We collect data every 10 minutes over the course of 240 min.

**Linear microrheology.** We determine linear viscoelastic properties from the thermal fluctuations of trapped microspheres, measured by recording the associated laser deflections for 100 seconds (Fig. 1B and D). We extract the elastic modulus  $G'(\omega)$  and viscous modulus  $G''(\omega)$ , from the thermal fluctuations using the generalized Stokes–Einstein relation (GSER) as described in ref. 81–83. In brief, we compute the normalized mean-squared displacements ( $\pi(\tau) = r^2(\tau)/2r^2$ ) of the thermal forces, which we convert into the Fourier domain *via*:

$$-\omega^2 \pi(\omega) = (1 - e^{-i\omega\tau_1}) \frac{\pi(\tau_1)}{\tau_1} + \dot{\pi}_\infty e^{-i\omega\tau_N} + \sum_{k=2}^N \left( \frac{\tau_k - \tau_{k-1}}{\tau_k - \tau_{k-1}} \right) (e^{-i\omega\tau_{k-1}} - e^{-i\omega\tau_k}),$$

where  $\tau$ , 1 and  $N$  represent the lag time and the first and last point of the oversampled  $\pi(\tau)$ .  $\dot{\pi}_\infty$  is the extrapolated slope of  $\pi(\tau)$  at infinity. Oversampling is done using the MATLAB function PCHIP.  $\pi(\omega)$  is related to viscoelastic moduli *via*:

$$G^*(\omega) = G'(\omega) + iG''(\omega) = \left( \frac{k}{6\pi R} \right) \left( \frac{1}{i\omega\pi(\omega)} - 1 \right),$$

where  $R$  and  $k$  represent the microsphere radius and trap stiffness. From  $G'(\omega)$  and  $G''(\omega)$  (Fig. S3, ESI†), we compute the magnitude of the complex modulus  $|G^*(\omega)| = [(G'(\omega))^2 + (G''(\omega))^2]^{1/2}$ , as well as the corresponding fractional elastic modulus  $g'(\omega) = G'(\omega)/|G^*(\omega)|$  (Fig. 2A and B) and complex viscosity  $\eta^*(\omega) = |G^*(\omega)|/\omega$  (Fig. 2C and D).

**Nonlinear microrheology.** We perform nonlinear microrheology measurements by displacing a trapped microsphere embedded in the sample through a distance  $s = 30 \mu\text{m}$  at a speed  $v = 60 \mu\text{m s}^{-1}$  using a piezoelectric nanopositioning stage (Mad City Laboratories) to move the sample relative to the microsphere (Fig. 1B and E), after which we hold the bead fixed for  $15 \text{ s}$  (Fig. 1C and F). We measure the force exerted on the microsphere for the duration of the measurement to determine

the force  $F(x, t)$  the solution exerts to resist the perturbation (Fig. 1E), and the relaxation of that force  $F_R(t)$  as the system relaxes to a new mechanically-steady-state (Fig. 1F).

### Differential dynamic microscopy (DDM)

We image fluorescent-labeled DNA molecules embedded in each topologically-active DNA solution with an Olympus IX73 epifluorescence microscope with an X-Cite LED light source, 488 nm emission and 535 nm detection filter cube, and a 60 $\times$ , 1.2NA objective (Olympus). Every 10 minutes after adding the restriction endonuclease ( $t_a = 0$ ), we capture times-series of images of 1440  $\times$  1920 pixels, providing a field-of-view of 87.3  $\mu\text{m} \times$  116  $\mu\text{m}$ . We use a Hamamatsu Orca-Flash 2.8 to capture images at a frame rate of 20 fps for 80 s. We divide each time-series into 25 separate 256  $\times$  256 pixel regions of interest (ROIs).

We perform DDM on each ROI using the PyDDM package on Github.<sup>72</sup> Following the method of DDM originally described by Cerbino and Trappe,<sup>70</sup> we compute the Fourier transform of the differences between images separated by a given lag time,  $\Delta t$ , for lag times that range from 50 ms (the time between frames) to 10 s. From this DDM analysis, we obtain the image structure function, which we model as  $D(q, \Delta t) = A(q)[1 - f(q, \Delta t)] + B$  where  $A(q)$  is the amplitude,  $B$  is the background, and  $f(q, \Delta t)$  is the intermediate scattering function (ISF). The amplitude and background terms depend on the properties of the imaging system, structure of the sample being imaged, and the camera noise. The ISF contains information about the system dynamics. The ISFs are fit to an exponential function of the form  $f(q, \Delta t) = \exp(-(\Delta t/\tau(q))^{s(q)})$  where  $\tau(q)$  is the characteristic decay time and  $s(q)$  is a stretching exponent that remained near or just under 1 for our data. For each ROI of the videos recorded every 10 minutes from  $t_a = 10$  min to  $t_a = 240$  min, we determine  $\tau(q)$ . For each  $t_a$ , we average  $\tau(q)$  across the 25 ROIs to determine the diffusion coefficient according to the equation  $D = 1/\tau q^2$ . We compute this diffusion coefficient for  $q$  values from 1 to 2.5  $\mu\text{m}^{-1}$ .

### Single-molecule tracking

To track the diffusion of single DNA molecules within topologically-active solutions, we add labeled DNA at a 10 $\times$  lower concentration compared to samples prepared for DDM analysis. The remaining imaging details and specifications are identical as described above. In every frame  $\sim$ 200 individual DNA molecules are visible and sufficiently separated and dispersed throughout the field-of-view to allow for accurate single-molecule tracking. We use custom particle tracking scripts (Python) to track the center-of-mass of individual DNA molecules and measure their  $x$  and  $y$  displacements ( $\Delta x$ ,  $\Delta y$ ) for a given lagtime  $\Delta t$ .<sup>84,85</sup> We then compute mean-squared displacements  $(\Delta x(\Delta t))^2$  and  $(\Delta y(\Delta t))^2$  from which we compute an average mean-squared displacement  $\text{MSD} = \frac{1}{2}((\Delta x)^2 + (\Delta y)^2)$  as a function of lag time. For MSDs that scale linearly with lag time, *i.e.*,  $\text{MSD} \sim (\Delta t)^\alpha$ , where  $\alpha = 1$ , the diffusion coefficient can be computed *via* the relation  $\text{MSD} = 2D\Delta t$ . Using this

approach, we determine DNA diffusion coefficients by fitting our measured MSDs from  $\Delta t = 0.1$  s to  $\Delta t = 0.7$  s, where the data is reliably linear, to a linear function and equating the slope to  $2D$ .

### Time-resolved Gel electrophoresis

To characterize the topological state of DNA solutions during linearization and fragmentation of circular DNA *via* HaeIII, we use direct current agarose gel electrophoresis to separate the different topologies and lengths of DNA. Specifically, we prepare 40  $\mu\text{L}$  samples as described above and incubate at RT for 240 min, during which we remove a 1  $\mu\text{L}$  aliquot from the reaction every 10 minutes and quench it with TE (100 mM Tris-HCl (pH 8), 1 mM EDTA) buffer and gel loading dye. We load 50 ng of DNA from each 'kinetic aliquot' onto a 1% agarose gel prepared with TAE (Tris-Acetate-EDTA) buffer, as described previously.<sup>14,65</sup> We run the gel at 5 V  $\text{cm}^{-1}$  for 2.5 hours, allowing for separation of the DNA into distinct bands corresponding to ring, supercoiled, full-length linear and fragmented linear (Fig. 6A). We use the standard  $\lambda$ -HindIII molecular marker to calibrate the gel to determine the DNA topology and length corresponding to each distinct band. To quantify the fraction of each topology at each 10 min interval, we perform band intensity analysis and use intensity as a proxy for mass. Specifically, we use a Life Technologies E-Gel Imager to image DNA bands on the gels and use Gel Quant Express software to perform image intensity analysis to determine the relative DNA mass in each band of a given lane.<sup>13</sup> Because HindIII produces so many fragments that are close in length, they cannot be well-separated or resolved, so we are unable to accurately determine the fraction of each fragment length. Instead, we treat all fragments as a single entity and compute the total mass fraction of all fragments together and compare to the mass fraction of full-length linear, ring and supercoiled DNA in each lane (Fig. 6A and B).

### Author contributions

R. M. R. A. designed the research and supervised experiments. P. N. and R. M. R. A. wrote the paper and interpreted data. P. N. and N. C. conducted experiments and analyzed data. R. M. supervised experiments, and helped with data analysis, interpretation, and writing the paper.

### Conflicts of interest

There are no conflicts to declare.

### Acknowledgements

P. N. and R. M. acknowledge the Research Corporation for Scientific Advancement for support. R. M. R. A. acknowledges the US Air Force Office of Scientific Research (Grant no. AFOSR-FA9550-17-1-0249) for support.

## References

- 1 B. Alberts, *Molecular Biology of the Cell*, Garland, 2004.
- 2 Q. Wen and P. A. Janmey, *Polymer Physics of the Cytoskeleton*, *Curr. Opin. Solid State Mater. Sci.*, 2011, **15**(5), 177–182, DOI: [10.1016/j.cossms.2011.05.002](https://doi.org/10.1016/j.cossms.2011.05.002).
- 3 R. L. P. Adams, J. T. Knowler and D. P. Leader, *The Biochemistry of the Nucleic Acids*, Springer Science & Business Media, 2013.
- 4 K. Liu, A. E. Patteson, E. J. Banigan and J. M. Schwarz, Dynamic Nuclear Structure Emerges from Chromatin Cross-Links and Motors, *Phys. Rev. Lett.*, 2021, **126**(15), 158101, DOI: [10.1103/PhysRevLett.126.158101](https://doi.org/10.1103/PhysRevLett.126.158101).
- 5 J. R. C. van der Maarel, *Introduction to Biopolymer Physics*, World Scientific, 2008.
- 6 A. F. Pegoraro, P. Janmey and D. A. Weitz, Mechanical Properties of the Cytoskeleton and Cells, *Cold Spring Harbor Perspect. Biol.*, 2017, **9**(11), a022038, DOI: [10.1101/cshperspect.a022038](https://doi.org/10.1101/cshperspect.a022038).
- 7 D. A. Fletcher and R. D. Mullins, Cell Mechanics and the Cytoskeleton, *Nature*, 2010, **463**(7280), 485–492, DOI: [10.1038/nature08908](https://doi.org/10.1038/nature08908).
- 8 S. Banerjee, M. L. Gardel and U. S. Schwarz, The Actin Cytoskeleton as an Active Adaptive Material, *Annu. Rev. Condens. Matter Phys.*, 2020, **11**(1), 421–439, DOI: [10.1146/annurev-conmatphys-031218-013231](https://doi.org/10.1146/annurev-conmatphys-031218-013231).
- 9 P. Chandrakar, J. Berezney, B. Lemma, B. Hishamunda, A. Berry, K.-T. Wu, R. Subramanian, J. Chung, D. Needleman, J. Gelles and Z. Dogic, Engineering Stability, Longevity, and Miscibility of Microtubule-Based Active Fluids, *Soft Matter*, 2022, **18**(9), 1825–1835, DOI: [10.1039/D1SM01289D](https://doi.org/10.1039/D1SM01289D).
- 10 R. J. McGorty, C. J. Currie, J. Michel, M. Sasanpour, C. Gunter, K. A. Lindsay, M. J. Rust, P. Katira, M. Das, J. L. Ross and R. M. Robertson-Anderson, Kinesin and Myosin Motors Compete to Drive Rich Multiphase Dynamics in Programmable Cytoskeletal Composites, *PNAS Nexus*, 2023, **2**(8), pgad245, DOI: [10.1093/pnasnexus/pgad245](https://doi.org/10.1093/pnasnexus/pgad245).
- 11 G. Lee, G. Leech, M. J. Rust, M. Das, R. J. McGorty, J. L. Ross and R. M. Robertson-Anderson, Myosin-Driven Actin-Microtubule Networks Exhibit Self-Organized Contractile Dynamics, *Sci. Adv.*, 2021, **7**(6), eabe4334, DOI: [10.1126/sciadv.abe4334](https://doi.org/10.1126/sciadv.abe4334).
- 12 B. A. Krajina, A. Zhu, S. C. Heilshorn and A. J. Spakowitz, Active DNA Olympic Hydrogels Driven by Topoisomerase Activity, *Phys. Rev. Lett.*, 2018, **121**(14), 148001, DOI: [10.1103/PhysRevLett.121.148001](https://doi.org/10.1103/PhysRevLett.121.148001).
- 13 D. Michieletto, P. Neill, S. Weir, D. Evans, N. Crist, V. A. Martinez and R. M. Robertson-Anderson, Topological Digestion Drives Time-Varying Rheology of Entangled DNA Fluids, *Nat. Commun.*, 2022, **13**(1), 4389, DOI: [10.1038/s41467-022-31828-w](https://doi.org/10.1038/s41467-022-31828-w).
- 14 J. Marfai, R. J. McGorty and R. M. Robertson-Anderson, Cooperative Rheological State-Switching of Enzymatically-Driven Composites of Circular DNA And Dextran, *Adv. Mater.*, 2023, 2305824, DOI: [10.1002/adma.202305824](https://doi.org/10.1002/adma.202305824).
- 15 A. Bonato, D. Marenduzzo, D. Michieletto and E. Orlandini, Topological Gelation of Reconnecting Polymers, *Proc. Natl. Acad. Sci. U. S. A.*, 2022, **119**(44), e2207728119, DOI: [10.1073/pnas.2207728119](https://doi.org/10.1073/pnas.2207728119).
- 16 D. Michieletto, Y. A. G. Fosado, E. Melas, M. Baiesi, L. Tubiana and E. Orlandini, Dynamic and Facilitated Binding of Topoisomerase Accelerates Topological Relaxation, *Nucleic Acids Res.*, 2022, **50**(8), 4659–4668, DOI: [10.1093/nar/gkac260](https://doi.org/10.1093/nar/gkac260).
- 17 Y. S. Kim, B. Kundukad, A. Allahverdi, L. Nordensköld, P. S. Doyle and J. R. C. Van Der Maarel, Gelation of the Genome by Topoisomerase II Targeting Anticancer Agents, *Soft Matter*, 2013, **9**(5), 1656–1663, DOI: [10.1039/C2SM27229F](https://doi.org/10.1039/C2SM27229F).
- 18 B. Kundukad and J. R. C. Van Der Maarel, Control of the Flow Properties of DNA by Topoisomerase II and Its Targeting Inhibitor, *Biophys. J.*, 2010, **99**(6), 1906–1915, DOI: [10.1016/j.bpj.2010.07.013](https://doi.org/10.1016/j.bpj.2010.07.013).
- 19 A. Pingoud, G. G. Wilson and W. Wende, Type II Restriction Endonucleases—A Historical Perspective and More, *Nucleic Acids Res.*, 2014, **42**(12), 7489–7527, DOI: [10.1093/nar/gku447](https://doi.org/10.1093/nar/gku447).
- 20 M. Gellert, DNA Topoisomerases, *Annu. Rev. Biochem.*, 1981, **50**(1), 879–910, DOI: [10.1146/annurev.bi.50.070181.004311](https://doi.org/10.1146/annurev.bi.50.070181.004311).
- 21 R. M. Robertson, S. Laib and D. E. Smith, Diffusion of Isolated DNA Molecules: Dependence on Length and Topology, *Proc. Natl. Acad. Sci. U. S. A.*, 2006, **103**(19), 7310–7314, DOI: [10.1073/pnas.0601903103](https://doi.org/10.1073/pnas.0601903103).
- 22 M. Rubinstein and R. H. Colby, *Polymer Physics*, Oxford University Press, Oxford, New York, 2003.
- 23 M. Doi and S. F. Edwards, *The Theory of Polymer Dynamics*, Oxford University Press, Oxford, 1986.
- 24 D. Michieletto, N. Nahali and A. Rosa, Glassiness and Heterogeneous Dynamics in Dense Solutions of Ring Polymers, *Phys. Rev. Lett.*, 2017, **119**(19), 197801, DOI: [10.1103/PhysRevLett.119.197801](https://doi.org/10.1103/PhysRevLett.119.197801).
- 25 P. Chaudhuri, L. Berthier and W. Kob, Universal Nature of Particle Displacements Close to Glass and Jamming Transitions, *Phys. Rev. Lett.*, 2007, **99**(6), 060604, DOI: [10.1103/PhysRevLett.99.060604](https://doi.org/10.1103/PhysRevLett.99.060604).
- 26 J. Wang and T. Ge, Crazing Reveals an Entanglement Network in Glassy Ring Polymers, *Macromolecules*, 2021, **54**(16), 7500–7511, DOI: [10.1021/acs.macromol.1c01080](https://doi.org/10.1021/acs.macromol.1c01080).
- 27 T. Ge, S. Panyukov and M. Rubinstein, Self-Similar Conformations and Dynamics in Entangled Melts and Solutions of Nonconcatenated Ring Polymers, *Macromolecules*, 2016, **49**(2), 708–722, DOI: [10.1021/acs.macromol.5b02319](https://doi.org/10.1021/acs.macromol.5b02319).
- 28 D. Parisi, S. Costanzo, Y. Jeong, J. Ahn, T. Chang, D. Vlassopoulos, J. D. Halverson, K. Kremer, T. Ge, M. Rubinstein, G. S. Grest, W. Srinin and A. Y. Grosberg, Nonlinear Shear Rheology of Entangled Polymer Rings, *Macromolecules*, 2021, **54**(6), 2811–2827, DOI: [10.1021/acs.macromol.0c02839](https://doi.org/10.1021/acs.macromol.0c02839).
- 29 Q. Huang, J. Ahn, D. Parisi, T. Chang, O. Hassager, S. Panyukov, M. Rubinstein and D. Vlassopoulos, Unexpected Stretching of Entangled Ring Macromolecules, *Phys. Rev. Lett.*, 2019, **122**(20), 208001, DOI: [10.1103/PhysRevLett.122.208001](https://doi.org/10.1103/PhysRevLett.122.208001).
- 30 D. Michieletto and T. Sakaue, Dynamical Entanglement and Cooperative Dynamics in Entangled Solutions of Ring and

- Linear Polymers, *ACS Macro Lett.*, 2021, **10**(1), 129–134, DOI: [10.1021/acsmacrolett.0c00551](https://doi.org/10.1021/acsmacrolett.0c00551).
- 31 M. Kapnistos, M. Lang, D. Vlassopoulos, W. Pyckhout-Hintzen, D. Richter, D. Cho, T. Chang and M. Rubinstein, Unexpected Power-Law Stress Relaxation of Entangled Ring Polymers, *Nat. Mater.*, 2008, **7**(12), 997–1002.
- 32 I. Chubak, E. Locatelli and C. N. Likos, Ring Polymers Are Much Stronger Depleting Agents than Linear Ones, *Mol. Phys.*, 2018, **116**(21–22), 2911–2926, DOI: [10.1080/00268976.2018.1503744](https://doi.org/10.1080/00268976.2018.1503744).
- 33 E. Ghobadpour, M. Kolb, M. R. Ejtehadi and R. Everaers, Monte Carlo Simulation of a Lattice Model for the Dynamics of Randomly Branching Double-Folded Ring Polymers, *Phys. Rev. E*, 2021, **104**(1), 014501, DOI: [10.1103/PhysRevE.104.014501](https://doi.org/10.1103/PhysRevE.104.014501).
- 34 J. Smrek, K. Kremer and A. Rosa, Threading of Unconcatenated Ring Polymers at High Concentrations: Double-Folded vs Time-Equilibrated Structures, *ACS Macro Lett.*, 2019, **8**(2), 155–160, DOI: [10.1021/acsmacrolett.8b00828](https://doi.org/10.1021/acsmacrolett.8b00828).
- 35 J. Smrek, J. Garamella, R. Robertson-Anderson and D. Michieletto, Topological Tuning of DNA Mobility in Entangled Solutions of Supercoiled Plasmids, *Sci. Adv.*, 2021, **7**(20), eabf9260, DOI: [10.1126/sciadv.abf9260](https://doi.org/10.1126/sciadv.abf9260).
- 36 M. A. Ubertini, J. Smrek and A. Rosa, Entanglement Length Scale Separates Threading from Branching of Unknotted and Non-Concatenated Ring Polymers in Melts, *Macromolecules*, 2022, **55**(23), 10723–10736, DOI: [10.1021/acs.macromol.2c01264](https://doi.org/10.1021/acs.macromol.2c01264).
- 37 C. Micheletti, I. Chubak, E. Orlandini and J. Smrek, Topology-Based Detection and Tracking of Deadlocks Reveal Aging of Active Ring Melts, *ACS Macro Lett.*, 2024, 124–129, DOI: [10.1021/acsmacrolett.3c00567](https://doi.org/10.1021/acsmacrolett.3c00567).
- 38 D. Michieletto, On the Tree-like Structure of Rings in Dense Solutions, *Soft Matter*, 2016, **12**(47), 9485–9500, DOI: [10.1039/C6SM02168A](https://doi.org/10.1039/C6SM02168A).
- 39 T. C. O'Connor, T. Ge, M. Rubinstein and G. S. Grest, Topological Linking Drives Anomalous Thickening of Ring Polymers in Weak Extensional Flows, *Phys. Rev. Lett.*, 2020, **124**(2), 027801, DOI: [10.1103/PhysRevLett.124.027801](https://doi.org/10.1103/PhysRevLett.124.027801).
- 40 G. S. Grest, T. Ge, S. J. Plimpton, M. Rubinstein and T. C. O'Connor, Entropic Mixing of Ring/Linear Polymer Blends, *ACS Polym. Au*, 2023, **3**(2), 209–216, DOI: [10.1021/acspolymersau.2c00050](https://doi.org/10.1021/acspolymersau.2c00050).
- 41 J. D. Halverson, W. B. Lee, G. S. Grest, A. Y. Grosberg and K. Kremer, Molecular Dynamics Simulation Study of Nonconcatenated Ring Polymers in a Melt. II. Dynamics, *J. Chem. Phys.*, 2011, **134**(20), 204905, DOI: [10.1063/1.3587138](https://doi.org/10.1063/1.3587138).
- 42 A. Y. Grosberg, Annealed Lattice Animal Model and Flory Theory for the Melt of Non-Concatenated Rings: Towards the Physics of Crumpling, *Soft Matter*, 2013, **10**(4), 560–565, DOI: [10.1039/C3SM52805G](https://doi.org/10.1039/C3SM52805G).
- 43 B. V. S. Iyer, S. Shanbhag, V. A. Juvekar and A. K. Lele, Self-Diffusion Coefficient of Ring Polymers in Semidilute Solution, *J. Polym. Sci., Part B: Polym. Phys.*, 2008, **46**(21), 2370–2379, DOI: [10.1002/polb.21569](https://doi.org/10.1002/polb.21569).
- 44 B. V. S. Iyer, A. K. Lele and V. A. Juvekar, Flexible Ring Polymers in an Obstacle Environment: Molecular Theory of Linear Viscoelasticity, *Phys. Rev. E: Stat., Nonlinear, Soft Matter Phys.*, 2006, **74**(2), 021805.
- 45 J. Smrek and A. Y. Grosberg, Minimal Surfaces on Unconcatenated Polymer Rings in Melt, *ACS Macro Lett.*, 2016, **5**(6), 750–754, DOI: [10.1021/acsmacrolett.6b00289](https://doi.org/10.1021/acsmacrolett.6b00289).
- 46 D. Parisi, M. Kaliva, S. Costanzo, Q. Huang, P. J. Lutz, J. Ahn, T. Chang, M. Rubinstein and D. Vlassopoulos, Nonlinear Rheometry of Entangled Polymeric Rings and Ring-Linear Blends, *J. Rheol.*, 2021, **65**(4), 695–711, DOI: [10.1122/8.0000186](https://doi.org/10.1122/8.0000186).
- 47 D. Kong, S. Banik, M. J. San Francisco, M. Lee, R. M. Robertson-Anderson, C. M. Schroeder and G. B. McKenna, Rheology of Entangled Solutions of Ring-Linear DNA Blends, *Macromolecules*, 2022, **55**(4), 1205–1217, DOI: [10.1021/acs.macromol.1c01672](https://doi.org/10.1021/acs.macromol.1c01672).
- 48 K. R. Peddireddy, M. Lee, C. M. Schroeder and R. M. Robertson-Anderson, Viscoelastic Properties of Ring-Linear DNA Blends Exhibit Nonmonotonic Dependence on Blend Composition, *Phys. Rev. Res.*, 2020, **2**(2), 023213, DOI: [10.1103/PhysRevResearch.2.023213](https://doi.org/10.1103/PhysRevResearch.2.023213).
- 49 D. Michieletto, D. Marenduzzo, E. Orlandini, G. P. Alexander and M. S. Turner, Threading Dynamics of Ring Polymers in a Gel, *ACS Macro Lett.*, 2014, **3**(3), 255–259, DOI: [10.1021/mz500060c](https://doi.org/10.1021/mz500060c).
- 50 D. G. Tsalirikis and V. G. Mavrantzas, Size and Diffusivity of Polymer Rings in Linear Polymer Matrices: The Key Role of Threading Events, *Macromolecules*, 2020, **53**(3), 803–820, DOI: [10.1021/acs.macromol.9b02099](https://doi.org/10.1021/acs.macromol.9b02099).
- 51 D. G. Tsalirikis, V. G. Mavrantzas and D. Vlassopoulos, Analysis of Slow Modes in Ring Polymers: Threading of Rings Controls Long-Time Relaxation, *ACS Macro Lett.*, 2016, **5**(6), 755–760, DOI: [10.1021/acsmacrolett.6b00259](https://doi.org/10.1021/acsmacrolett.6b00259).
- 52 K. R. Peddireddy, R. Clairmont and R. M. Robertson-Anderson, Polymer Threadings and Rigidity Dictate the Viscoelasticity of Entangled Ring-Linear Blends and Their Composites with Rigid Rod Microtubules, *J. Rheol.*, 2023, **67**(1), 125–138, DOI: [10.1122/8.0000529](https://doi.org/10.1122/8.0000529).
- 53 B. V. S. Iyer, A. K. Lele and S. Shanbhag, What Is the Size of a Ring Polymer in a Ring-Linear Blend?, *Macromolecules*, 2007, **40**(16), 5995–6000.
- 54 J. D. Halverson, G. S. Grest, A. Y. Grosberg and K. Kremer, Rheology of Ring Polymer Melts: From Linear Contaminants to Ring-Linear Blends, *Phys. Rev. Lett.*, 2012, **108**(3), 038301, DOI: [10.1103/PhysRevLett.108.038301](https://doi.org/10.1103/PhysRevLett.108.038301).
- 55 T. C. O'Connor, T. Ge and G. S. Grest, Composite Entanglement Topology and Extensional Rheology of Symmetric Ring-Linear Polymer Blends, *J. Rheol.*, 2022, **66**(1), 49–65, DOI: [10.1122/8.0000319](https://doi.org/10.1122/8.0000319).
- 56 D. Michieletto, D. Marenduzzo, E. Orlandini and M. S. Turner, Ring Polymers: Threadings, Knot Electrophoresis and Topological Glasses, *Polymers*, 2017, **9**(8), 349, DOI: [10.3390/polym9080349](https://doi.org/10.3390/polym9080349).
- 57 K. R. Peddireddy, M. Lee, Y. Zhou, S. Adalbert, S. Anderson, C. M. Schroeder and R. M. Robertson-Anderson,

- Unexpected Entanglement Dynamics in Semidilute Blends of Supercoiled and Ring DNA, *Soft Matter*, 2020, **16**(1), 152–161, DOI: [10.1039/C9SM01767D](https://doi.org/10.1039/C9SM01767D).
- 58 S. Banik, D. Kong, M. J. San Francisco and G. B. McKenna, Monodisperse Lambda DNA as a Model to Conventional Polymers: A Concentration-Dependent Scaling of the Rheological Properties, *Macromolecules*, 2021, **54**(18), 8632–8654, DOI: [10.1021/acs.macromol.0c02537](https://doi.org/10.1021/acs.macromol.0c02537).
- 59 R. M. Robertson and D. E. Smith, Self-Diffusion of Entangled Linear and Circular DNA Molecules: Dependence on Length and Concentration, *Macromolecules*, 2007, **40**(9), 3373–3377.
- 60 T. M. Squires and T. G. Mason, Fluid Mechanics of Microrheology, *Annu. Rev. Fluid Mech.*, 2010, **42**(1), 413–438, DOI: [10.1146/annurev-fluid-121108-145608](https://doi.org/10.1146/annurev-fluid-121108-145608).
- 61 R. M. Robertson and D. E. Smith, Direct Measurement of the Confining Forces Imposed on a Single Molecule in a Concentrated Solution of Circular Polymers, *Macromolecules*, 2007, **40**(24), 8737–8741, DOI: [10.1021/ma071440e](https://doi.org/10.1021/ma071440e).
- 62 S.-Q. Wang, S. Ravindranath, Y. Wang and P. Boukany, New Theoretical Considerations in Polymer Rheology: Elastic Breakdown of Chain Entanglement Network, *J. Chem. Phys.*, 2007, **127**(6), 064903, DOI: [10.1063/1.2753156](https://doi.org/10.1063/1.2753156).
- 63 C. D. Chapman and R. M. Robertson-Anderson, Nonlinear Microrheology Reveals Entanglement-Driven Molecular-Level Viscoelasticity of Concentrated DNA, *Phys. Rev. Lett.*, 2014, **113**(9), 098303, DOI: [10.1103/PhysRevLett.113.098303](https://doi.org/10.1103/PhysRevLett.113.098303).
- 64 M. Khan, K. Regan and R. M. Robertson-Anderson, Optical Tweezers Microrheology Maps the Dynamics of Strain-Induced Local Inhomogeneities in Entangled Polymers, *Phys. Rev. Lett.*, 2019, **123**(3), 038001, DOI: [10.1103/PhysRevLett.123.038001](https://doi.org/10.1103/PhysRevLett.123.038001).
- 65 D. Michieletto, P. Neill, S. Weir, D. Evans, N. Crist, V. A. Martinez and R. M. Robertson-Anderson, Topological Digestion Drives Time-Varying Rheology of Entangled DNA Fluids, *Nat. Commun.*, 2022, **13**(1), 4389, DOI: [10.1038/s41467-022-31828-w](https://doi.org/10.1038/s41467-022-31828-w).
- 66 R. M. Robertson and D. E. Smith, Direct Measurement of the Intermolecular Forces Confining a Single Molecule in an Entangled Polymer Solution, *Phys. Rev. Lett.*, 2007, **99**(12), 126001, DOI: [10.1103/PhysRevLett.99.126001](https://doi.org/10.1103/PhysRevLett.99.126001).
- 67 M. Doi and S. F. Edwards, *The Theory of Polymer Dynamics*, Clarendon Press, 1988.
- 68 C. D. Chapman, K. Lee, D. Henze, D. E. Smith and R. M. Robertson-Anderson, Onset of Non-Continuum Effects in Microrheology of Entangled Polymer Solutions, *Macromolecules*, 2014, **47**(3), 1181–1186, DOI: [10.1021/ma401615m](https://doi.org/10.1021/ma401615m).
- 69 R. Cerbino, F. Giavazzi and M. E. Helgeson, Differential Dynamic Microscopy for the Characterization of Polymer Systems, *J. Polym. Sci.*, 2022, **60**(7), 1079–1089, DOI: [10.1002/pol.20210217](https://doi.org/10.1002/pol.20210217).
- 70 R. Cerbino and V. Trappe, Differential Dynamic Microscopy: Probing Wave Vector Dependent Dynamics with a Microscope, *Phys. Rev. Lett.*, 2008, **100**(18), 188102, DOI: [10.1103/PhysRevLett.100.188102](https://doi.org/10.1103/PhysRevLett.100.188102).
- 71 F. Giavazzi, D. Brogioli, V. Trappe, T. Bellini and R. Cerbino, Scattering Information Obtained by Optical Microscopy: Differential Dynamic Microscopy and Beyond, *Phys. Rev. E: Stat., Nonlinear, Soft Matter Phys.*, 2009, **80**(3), 031403, DOI: [10.1103/PhysRevE.80.031403](https://doi.org/10.1103/PhysRevE.80.031403).
- 72 H. N. Verwei, G. Lee, G. Leech, I. I. Petitjean, G. H. Koenderink, R. M. Robertson-Anderson and R. J. McGorty, Quantifying Cytoskeleton Dynamics Using Differential Dynamic Microscopy, *JoVE J. Vis. Exp.*, 2022, **184**, e63931, DOI: [10.3791/63931](https://doi.org/10.3791/63931).
- 73 J. Garamella, K. Regan, G. Aguirre, R. J. McGorty and R. M. Robertson-Anderson, Anomalous and Heterogeneous DNA Transport in Biomimetic Cytoskeleton Networks, *Soft Matter*, 2020, **16**(27), 6344–6353, DOI: [10.1039/d0sm00544d](https://doi.org/10.1039/d0sm00544d).
- 74 P. Khanal, K. R. Peddireddy, J. Marfai, R. McGorty and R. M. Robertson-Anderson, DNA Topology Dictates Emergent Bulk Elasticity and Hindered Macromolecular Diffusion in DNA-Dextran Composites, *J. Rheol.*, 2022, **66**(4), 699–715, DOI: [10.1122/8.0000447](https://doi.org/10.1122/8.0000447).
- 75 R. M. Robertson and D. E. Smith, Strong Effects of Molecular Topology on Diffusion of Entangled DNA Molecules, *Proc. Natl. Acad. Sci. U. S. A.*, 2007, **104**(12), 4824–4827, DOI: [10.1073/pnas.0700137104](https://doi.org/10.1073/pnas.0700137104).
- 76 R. M. Robertson, S. Laib and D. E. Smith, Diffusion of Isolated DNA Molecules: Dependence on Length and Topology, *Proc. Natl. Acad. Sci. U. S. A.*, 2006, **103**(19), 7310–7314, DOI: [10.1073/pnas.0601903103](https://doi.org/10.1073/pnas.0601903103).
- 77 S. Laib, R. M. Robertson and D. E. Smith, Preparation and Characterization of a Set of Linear DNA Molecules for Polymer Physics and Rheology Studies, *Macromolecules*, 2006, **39**(12), 4115–4119, DOI: [10.1021/ma060146a](https://doi.org/10.1021/ma060146a).
- 78 S. Pan, D. At Nguyen, T. Sridhar, P. Sunthar and J. Ravi Prakash, Universal Solvent Quality Crossover of the Zero Shear Rate Viscosity of Semidilute DNA Solutions, *J. Rheol.*, 2014, **58**(2), 339–368, DOI: [10.1122/1.4861072](https://doi.org/10.1122/1.4861072).
- 79 J. R. Prakash, Universal Dynamics of Dilute and Semidilute Solutions of Flexible Linear Polymers, *Curr. Opin. Colloid Interface Sci.*, 2019, **43**, 63–79, DOI: [10.1016/j.cocis.2019.03.001](https://doi.org/10.1016/j.cocis.2019.03.001).
- 80 R. M. Robertson-Anderson, Optical Tweezers Microrheology: From the Basics to Advanced Techniques and Applications, *ACS Macro Lett.*, 2018, **7**(8), 968–975, DOI: [10.1021/acsmacrolett.8b00498](https://doi.org/10.1021/acsmacrolett.8b00498).
- 81 M. Tassieri, R. M. L. Evans, R. L. Warren, N. J. Bailey and J. M. Cooper, Microrheology with Optical Tweezers: Data Analysis, *New J. Phys.*, 2012, **14**(11), 115032, DOI: [10.1088/1367-2630/14/11/115032](https://doi.org/10.1088/1367-2630/14/11/115032).
- 82 M. Tassieri, G. M. Gibson, R. M. L. Evans, A. M. Yao, R. Warren, M. J. Padgett and J. M. Cooper, Measuring Storage and Loss Moduli Using Optical Tweezers: Broadband Microrheology, *Phys. Rev. E: Stat., Nonlinear, Soft Matter Phys.*, 2010, **81**(2), 026308, DOI: [10.1103/PhysRevE.81.026308](https://doi.org/10.1103/PhysRevE.81.026308).
- 83 T. G. Mason, Estimating the Viscoelastic Moduli of Complex Fluids Using the Generalized Stokes–Einstein Equation | SpringerLink, *Rheol. Acta*, 2000, **39**, 371–378.



- 84 S. J. Anderson, C. Matsuda, J. Garamella, K. R. Peddireddy, R. M. Robertson-Anderson and R. McGorty, Filament Rigidity Vies with Mesh Size in Determining Anomalous Diffusion in Cytoskeleton, *Biomacromolecules*, 2019, **20**(12), 4380–4388, DOI: [10.1021/acs.biomac.9b01057](https://doi.org/10.1021/acs.biomac.9b01057).
- 85 J. Y. Sheung, J. Garamella, S. K. Kahl, B. Y. Lee, R. J. McGorty and R. M. Robertson-Anderson, Motor-Driven Advection Competes with Crowding to Drive Spatiotemporally Heterogeneous Transport in Cytoskeleton Composites, *Front. Phys.*, 2022, **10**, 1055441.